

Pen Needles



*Pen injectors compatible with the pen needle:
 NovoPen®5, Autopen®24, Gansulin Pen, Ypsopen®, UNIPEN®, XiuLin Pen,
 HumaPen®ErgoII, Byetta, Victoza Pen
 * All the brands are trademarks of their respective owners.
 Please read insulin pen instructions for use.

- High Flow Rate
- Silicone Coating
- Classic bevel
- Gamma Irradiation

- Sterilized by Gamma Irradiation
- Non-Toxic, Non-Pyrogenic
- Without Natural Rubber Latex
- Full range sizes, Private Label Available

Ordering information

Reorder #	Color	Needle gauge	Needle length	Standard packing
01-3204		32G needle	4 mm	100 pcs/box
01-3105		31G needle	5 mm	2,000 pcs/case
01-3106		31G needle	6 mm	

1. Needle gauge refers to the size of the needle. The higher the gauge, the thinner the needle;
2. Please contact our sales person if you have special requests for an unique size.